Süddeutsche

Hauser

Chirurgie-Mechanik

GB

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Instructions for Use Stainless steel instruments

Symbols on device and packaging

Caution, general warning sign
Caution, consult accompanying documents

Scope

These instructions for use refer to instruments, of which all surfaces are directly accessible and visible after self-explanatory disassembly (e.g. screwing, inserting, pulling), manufactured by Max Hauser Süddeutsche Chirurgiemechanik GmbH. The instrument does not include parts made from nickel silver or tin and is neither nickel- nor chrome-plated.

Intended use

These instructions for use describe an instrument for universal surgical use in all specialist areas.

Safe handling and preparation

- ▶ The device and its accessories may only be used and applied by medical practitioners with the requisite training, knowledge and experience as improper use of an instrument may result in serious injury.
- ▶ Carefully read and follow the instructions for use and keep them within reach.
- ▶ Use the device only in accordance with its intended use.
- ▶ After removing the transport packaging, clean the new device prior to its initial (manual or automatic) sterilization.
- ▶ Safely store new/unused devices in a dry and clean place.
- ▶ Prior to each use, inspect the device for loose, bent, broken, cracked, worn, or damaged parts.
- ▶ Do not use the device if it is damaged or defective. Immediately eliminate the device if damaged.
- ▶ Immediately replace damaged parts with original spare parts.
- ▶ Improper use of an instrument may result in serious injury.

Reprocessing of the instrument

General safety information

- National statutory regulations, national and international standards and directives, and local, clinical hygiene instructions for reprocessing are to be complied with.
- For patients with Creutzfeldt-Jakob disease (CJD), suspected CJD, or possible variants of CJD, the relevant national regulations concerning the reprocessing of devices are to be complied with.
- It should be noted that successful reprocessing of this medical device can only be guaranteed following prior validation of the reprocessing method. The operator/reprocessor is responsible for this.

General cleaning information

Tissue, body fluids or other surgical residue that is left to dry on surgical instruments can prevent effective cleaning and may lead to corrosion. For that reason, reprocessing should begin within six (6) hours after surgery. Do not use fixating detergents, hot water >45° C or fixating disinfectants containing active ingredients such as aldehyde or alcohol.

Carefully follow the manufacturers' instructions when preparing cleaning agents. An overdose of neutralizing agents or of basic cleaners may result in chemical damage and/or fading. This may render laser marking illegible and/or electronically unrecognizable, even with instruments made of stainless steel.

Likewise, chlorine and chloride-containing residues (e.g. surgical residues, drugs, saline solutions and similar substances in the water that is used for cleaning, disinfection and sterilization) may cause corrosion damage (pitting, stress corrosion) and thus result in damage to the device. These residues must be removed by thorough rinsing with demineralized water and subsequent drying of the device.

Do not overload the washer unit. Ensure that no areas are left unwashed. Perform additional drying, if necessary.

Only use processing chemicals that have been tested and approved (e.g. VAH, FDA, or CE mark) and that are compatible with the device's materials according to the chemical manufacturers' recommendations. The chemical manufacturers' specifications must be strictly observed. Failure to do so may result in the following problems:

- ▶ Optical deterioration of materials, e.g. fading or discoloration on the surface
- ▶ Material damage, e.g. corrosion, cracks, fracturing, premature aging
- ▶ Do not use metal cleaning brushes or surface-damaging abrasive agents that might cause corrosion
- ▶ Further detailed advice on hygienically safe and material-/value-preserving reprocessing methods is available at www.a-k-i.org Publications Red Brochure – Proper reprocessing of instruments.

Preparation at the place of use

- ▶ Rinse surfaces not available for visible inspection, preferably with demineralized water, using a disposable syringe or water pressure gun.
- ▶ To the extent possible, remove visible surgical residues with a damp, lint-free cloth.
- ▶ Place the device in a sealed transportation container and take it to the facility's cleaning and disinfection unit as soon as possible. Reprocessing must be initiated within two (2) hours after use.

Cleaning and disinfection

Device-specific safety information on the reprocessing method

Caution

Inappropriate cleaning/disinfecting agents and/or excessive temperatures may cause damage to the device!

- ▶ Use cleaning agents and disinfectants approved for stainless steel in accordance with the manufacturer's instructions.
- ▶ Observe the instructions regarding concentration, temperature and application time.

- ▶ Do not exceed the maximum permitted cleaning temperature of 90 °C.

Validated cleaning and disinfection procedure

| Validated procedure | Specific requirements | Reference |
|--|---|--|
| Manual pre-cleaning with subsequent automated alkaline cleaning and thermal disinfection | <ul style="list-style-type: none"> ▪ suitable cleaning brush and water pressure gun ▪ disposable 20 ml syringe ▪ Put jaw covers on the device ▪ Place the device in an adequate sterilization tray (ensure that no areas can be left unwashed). | Chapter "Automatic washer/disinfector" with manual pre-cleaning and sub-chapter: <ul style="list-style-type: none"> ▪ Chapter "Manual Pre-Cleaning" ▪ Chapter "Automatic alkaline cleaning and thermal disinfection" |

Information on automatic cleaning and disinfection with manual pre-cleaning

- Only automatic washers and disinfectors with certified effectiveness (e.g. FDA approval or CE mark according to DIN EN ISO 15883) may be used.
- Washers and disinfectors must be serviced and checked at regular intervals.

Manual Pre-Cleaning

| Phase | Step | D (°C/°F) | t (min) | Conc. (%) | Water quality | Chemistry |
|-------|-----------------------|-----------|---------|-----------|---------------|--|
| I | Disinfecting cleaning | RT (cold) | >15 | 2 | DW | Concentrate, aldehyde-free, phenol-free, pH ~ 9* |
| II | Rinsing | RT (cold) | 1 | - | DW | - |

DW: Drinking water

RT Room temperature

- ▶ Disassemble prior to cleaning

Phase I

- ▶ Start by removing coarse contaminations from the device using a suitable water pressure gun.
- ▶ Fully immerse the device in the cleaning/disinfecting solution for a minimum of 15 minutes. Ensure that all accessible surfaces are wet.
- ▶ Clean the device immersed in the solution with a suitable cleaning brush until any visible contamination has been removed.
- ▶ If applicable, clean non-visible surfaces for a minimum of 1 min using a suitable cleaning gun.
- ▶ During the cleaning process, activate non-rigid component connections such as adjustable screws, hinges, etc., by opening and closing or moving back and forth, as appropriate, and
- ▶ Subsequently rinse these connections thoroughly with the cleaning / disinfection solution using a suitable disposable syringe. Repeat at least five (5) times.

Phase II

- ▶ Thoroughly rinse/flush all accessible instrument surfaces with a water pressure gun.
- ▶ Activate non-rigid component connections such as adjustable screws, hinges, etc., by opening and closing or moving back and forth during the cleaning process.

Automatic alkaline cleaning and thermal disinfection

Device type: Single-chamber washer-disinfector without ultrasound

| Phase | Step | T (°C/°F) | t (min) | Water quality | Chemistry |
|-------|-----------|-----------|---------|---------------|---|
| I | Pre-rinse | <25/77 | 3 | DW | - |
| II | Cleaning | 55/131 | 10 | FDW | <ul style="list-style-type: none"> ▪ Concentrate: <ul style="list-style-type: none"> - pH neutral (pH ~ 13 for stainless-steel instruments) - <5% anionic surfactant ▪ ready-to-use solution 0.5% |

| | | | | | |
|------|---|---------|---|-----|--|
| | | | | | - pH neutral (pH ~ 11 for stainless-steel instruments) |
| III | Intermediate rinse | > 10/50 | 1 | FDW | - |
| IV | Thermal disinfection | 90/194 | 5 | FDW | |
| V | Drying | - | - | - | In accordance with the program for the cleaning and disinfecting machine |
| DW: | Drinking water | | | | |
| FDW: | Fully desalinated water (demineralized, low microbiological contamination: drinking water quality at least) | | | | |

- ▶ Following automatic washing/disinfection, check visible surfaces for residual contamination.

Inspection, maintenance and checks

- ▶ Allow the device to cool down to room temperature.
- ▶ After each washing, disinfection and drying cycle, inspect the device to ensure it is dry, clean, functioning, undamaged, e.g. check for insulation damage, corrosion, loose, bent, broken, cracked, worn and detached parts.
- ▶ Dry the device if wet or damp.
- ▶ Repeat washing and disinfection steps for contaminated devices.
- ▶ Verify the device function.
- ▶ Immediately eliminate damaged or non-functional devices and send them to Max Hauser Süddeutsche Chirurgiemechnik GmbH for evaluation, clarification of the root cause and repair.
- ▶ Verify compatibility with associated devices.

Packaging

- ▶ Place the device in its dedicated holder or in a suitable tray. Ensure that sharp edges are covered.
- ▶ Adequately wrap the trays in accordance with the sterilization procedure used.
- ▶ Ensure that the packaging sufficiently protects the device from recontamination during storage.

Steam sterilization

- ▶ Ensure that the sterilizing agent can come into contact with all external and internal surfaces (e.g., by opening any valves, stopcocks, etc.)
- ▶ Validated sterilization procedure
 - Steam sterilization (fractional vacuum process)
 - Steam sterilizer according to DIN EN 285, validated against DIN EN ISO 17665
 - Sterilization using the fractionated vacuum process at 134°C, holding time 5 min
- ▶ If several devices are sterilized at the same time in the same steam sterilizer: Ensure that the maximum permitted load according to the manufacturers' specifications is not exceeded.

Storage

- ▶ Store sterile devices in germ-proof packaging, protected from dust, in a dry, dark, temperature-controlled room.

Technical Service

Caution

Risk of injury and/or malfunction!

- ▶ Modification of the device is not permitted.

Modifications of medical devices may result in a loss of potential guarantee/warranty claims and forfeiture of applicable licenses.

Manufacturer

Max Hauser Süddeutsche Chirurgiemechnik GmbH
Föhrenstr. 33
D-78532 Tuttlingen / Germany
Phone: +49 7461 770 564 0
Fax: +49 7461 770 564 99
E-mail: info@max-hauser.com

Service address

Max Hauser Süddeutsche Chirurgiemechnik GmbH
Föhrenstr. 33
D-78532 Tuttlingen / Germany
Phone: +49 7461 770 564 0
Fax: +49 7461 770 564 99
E-mail: info@max-hauser.com

Disposal

- ▶ Dispose or recycle the device in accordance with the applicable national regulations!

Symbols on the label

Manufacturer

Reference number

Batch number

Date of manufacture

Non-sterile – sterilize prior to use

Read the instructions for use

The device complies with the requirements of EU Directive 93/42/EEC for Class 1 products

0297 The device complies with the requirements of EU Directive 93/42/EEC

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Ürün, sınıf I ürünler için AB yönetmeliği 93/42/EWG gereksinimleri yerine getirmektedir
Ürün AB yönetmeliği 93/42/EWG gereksinimleri yerine getirmektedir.